

Coactivation Patterns Pipeline

Brain Connectivity Lab
University of Electronic Science and Technology of China
Chengdu 611731, China
david_young1994@foxmail.cpm

If you used our CAP pipeline scripts, please cite our paper

‘Yang, H., Zhang, H., Di, X., Wang, S., Meng, C., Tian, L., & Biswal, B. (2021). Reproducible coactivation patterns of functional brain networks reveal the aberrant dynamic state transition in schizophrenia. *NeuroImage*, 118193.’

0) Prerequisites

This CAP pipeline is based on MATLAB, please add the folder ‘**CAP_functions**’ using ‘**add with subfolders**’. There are three folders within the ‘**CAP_functions**’.

› work (E:) › CoactivationPatterns › CAP_functions ›

Name	Date modified
✓ File folder (3)	
Functions	2021/10/11 16:40
Miscellaneous	2021/10/11 16:41
Plots	2021/10/9 16:47


- ‘**Functions**’ is the main folder of the pipeline. It includes the k-means of CAPs, generates the NIFTI files for the group-averaged CAPs, and calculates temporal dynamic measures including fraction of time, persistence, counts, transition probability, and so on.
- ‘**Plots**’ includes scripts to plot the results.
- ‘**Miscellaneous**’ includes some supplementary scripts and files used in the pipeline.

In addition, **SPM12** (<https://www.fil.ion.ucl.ac.uk/spm/software/spm12/>) and **DPARSF** (<http://www.rfmri.org/DPARSF>) are required to import NIFTI files. **BrainNet Viewer** (<https://www.nitrc.org/projects/bnv>) is used to generate the CAP spatial maps. Please download and add these toolboxes to your MATLAB before using our pipeline.

1) Data structure

The data of HC and Patient should be saved in separate folders, named **HC** and **PAT**. Within the group folder, each subject also has one folder to save the preprocessed images (either 4D or 3D NIFTI files are ok).


For example, ‘E:\Data\CAP_project\HC\HC001\4D_HC001.nii’



2) The CAP pipeline

a) Main function

The main function is **'CAP_Framework'**, and it will call the **'CAP_Framework_HC'** and **'CAP_Framework_PAT'**. To run the CAP pipeline, you only need to set some parameters within the **'CAP_Framework'**, and run the **'CAP_Framework'**.



If you only have one HC group, then annotate the **CAP – PAT** part. As for the patient group, we don't perform the k-means clustering, and instead, we directly match volumes from patients' data with the group-averaged CAPs obtained from the HC group. Therefore, if you only have one patient group and without the control group, you can still annotate the **CAP – PAT** part, and use **'CAP_Framework_HC'** to generate the CAPs of patients and calculate the temporal dynamic measures.

```
% CAP - PAT, CAP pipeline for the Patient group
CAP_PAT_path = [working_dir 'Results' filesep 'PAT' filesep]; % create the PAT folder to save their CAP results
mkdir(CAP_PAT_path)

subjfolder_PAT = dir([Data_Path.Sub 'PAT']); % the data folder of PAT

% run the CAP pipeline in the Patient group
CAP_Framework_PAT(CAP_PAT_path, subjfolder_PAT, Data_Path, K_Means_Para)
```

b) the settings for the CAP analysis

To run the CAP pipeline, you only need to set the filepath, working directory and k-means parameters. Then run the **‘CAP_Framework’**, and it will finish generating the CAPs and calculating their temporal dynamic measures.

```
%% Default settings
clear
clc

Data_Path.Sub = 'E:\Data\CAP_project\': % set the subject folder
Data_Path.ROI = 'E:\CoactivationPatterns\CAP_functions\Miscellaneous\yeo_400_3mm.nii': % set the ROI template filepath for BOLD signal extraction

working_dir = 'E:\Project\CAP\': % set the working directory, all results and figures will be saved in this folder

% set the parameters for kmeans
K_Means_Para.Range = [2:1:21]: % the range of K, e.g., from 2 to 21 with step length = 1
K_Means_Para.Num = length(K_Means_Para.Range);
K_Means_Para.DistanceFunc = 'correlation': % distance function, e.g., 'sqeuclidean' | 'cityblock' | 'cosine' | 'correlation' | 'hamming'
K_Means_Para.ReplicateNum = 100: % replicate number
```

c) the settings for plots

There are several setting files for plots.

'\CAP_functions\Miscellaneous\CAP_colormap.mat' is the colormap to plot the matrix, e.g., the between-CAP similarity matrix.

'\CAP_functions\Plots\CAP_BrainNet.mat' is the setting option of BrainNet Viewer to plot the CAP spatial map. You can change it based on your requirement.

You also need to change the **‘brain_path’** to the filepath of your BrainNet Viewer. You can find and use other underlay in SurfTemplate.

```
%% Plot CAP spatial patterns
brain_path = 'E:\CoactivationPatterns\CAP_functions\Miscellaneous\BrainNetViewer_20191031\Data\SurfTemplate\BrainMesh_ICBM152_smoothed.mv';
setting_path = 'E:\CoactivationPatterns\CAP_functions\Plots\CAP_BrainNet.mat': % set the BrainNet Viewer options

for K_i = 1 : K_Means_Para.Num
    K = K_Means_Para.Range(K_i);
    CAP_str = ['CAP' num2str(K)];
    vol_path = [CAP_HC_path CAP_str]: % volume path for the group averaged CAPs
    outpath = [plot_path CAP_str];

    BrainNet_CAP(brain_path, setting_path, vol_path, outpath, K)
end
close all
```


3) Results and figures

For example, the CAP results of $K = 4$ in the HC group, their CAP spatial maps would be saved in the folder, **‘\working directory\Results\HC\CAP4’**


▼ NII File (12)

- Group_CAP_1_mean.nii
- Group_CAP_1_std.nii
- Group_CAP_1_z.nii
- Group_CAP_2_mean.nii
- Group_CAP_2_std.nii
- Group_CAP_2_z.nii
- Group_CAP_3_mean.nii
- Group_CAP_3_std.nii
- Group_CAP_3_z.nii
- Group_CAP_4_mean.nii
- Group_CAP_4_std.nii
- Group_CAP_4_z.nii


And their temporal dynamic measures would be saved in the ‘CAP_matrices.mat’.

 CAP_matrices.mat


The Figures of the CAP maps would be saved in the folder,
‘\working directory\Plots\HC\CAP4\’




Group_CAP_1_z.tif



Group_CAP_2_z.tif



Group_CAP_3_z.tif



Group_CAP_4_z.tif